Sterilization Case



Nozzle Spray Adaptor for Slim (P200-EZ-SM) distantion of Spray Adaptor for Bone Saw (P200-EZ-BM) Sprav Adaptor

EZ Spray Adaptor

E-Type Spray

bracket for the cable on the side.

Irrigation Tube



Z299

iCart-L Mobile Trolley

- Dimensions :
- W360 x D420 x H800 (mm) • Weight : 10 kg



PD-IT2S

Surgical Irrigation Tube • 5 pcs/ set (Sterilized)

Spray Mist Absorber

A convenient way to avoid mist spray from escaping into the environment. Simply insert the handpiece into the Absorber Pack then apply the spray lubricant.



Maintenance

PANA SPRAY Plus is a specially formulated, highly effective general handpiece lubricant. Frequent application of PANA SPRAY Plus will assure handpiece performance and durability. Apply PANA SPRAY Plus to handpieces after "Cleaning" in the maintenance process.

. Each handpiece

- Mailer

(P200-EZ-WPD





Z182100

PANA Spray Plus

- For Handpieces, 480 mL
- Alcohol-contained Lubricant

Y900630

- PANA Spray Plus (6pcs)
- Pack of 6



Z019090

E-Type Spray Nozzle

• For Handpieces

PD-EZ-M

EZ Spray Adaptor • PANA Spray Plus is not included





Spray Mist Absorber Pack of 10